

## **PYCNOMETER TOP AND GLASS JAR 48-D0441/A**

STANDARD: ASTM C128

Pycnometer top and 1000ml glass jar set for determining specific gravity of fine aggregate. Top is spun brass with 10mm hole in one end; threaded end fits 1000ml glass jar. Includes rubber gasket that fits on jar mouth to prevent fine particles from becoming deposited in the threads.



## **SAND ABSORPTION CONE AND TAMPER 48-D0440**

STANDARD: EN 12274-3

Cone and tamper manufactured according to the specifications and is also used in determining the specific gravity and absorption of fine aggregates.

SA-C Cone weight: 250 g  
SA-T Tamper weight: 430g

